

No image found for product type unknown

2 Pipette Stand (Verticle) –

Sku : GE6415



www.geomatrixlabs.com

Description

Geomatrix Labs stable and non-corrosive detachable stand capable of holding 28 pipettes with a maximum diameter of 14 mm. The heavy base measures 225 mm in diameter, and a Polypropylene-coated square rod measuring 375 mm in length is screwed into the base. The top discs, fixable on the rod, feature holes for holding pipettes, while the lower one has pilot holes for pipette tips and liquid drainage to the base. These discs can be adjusted on the rod based on pipette length. Ensure the arrows provided on the discs are kept parallel to center corresponding holes accurately.