

No image found for product type unknown

## Autoclave Vertical Double Walled – Radial Locking Modular Cabinet Type

**SKU :** GI1016



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Description

In steam pressure sterilization autoclaves are ideal for all applications requiring, assured and reliable total destruction of all living micro-organisms

Our top loading Autoclaves are the ideal companion for the modern microbiologists and institutes, pharma companies, hospitals etc.

### FEATURES :

- In Outer Rectangular body, inner, outer chamber and support ring are made of Stainless Steel (SS-304 grade)
- Heaters cover stand and carrier/bucket are made of stainless steel.
- Lid made of thick machined stainless steel plate.
- Radial locking lid ensures user safety and safe locking.
- Lid is fitted with pressure gauge, spring loaded safety valve, double safety valve, manual exhaust valve, vacuum breaker as a safety device and water level indicator.
- Easy portability & placement, unit is mounted on castor wheels for easy movement
- Drain valve facilitates easy cleaning of the chamber.
- Moulded joint less silicon rubber gasket.
- All internal joints are argon arc welded, ground and polished to give crevice free internals.
- Miniature circuit breaker for mains on/off plus safety of control panel.
- Industrial grade energy efficient heaters reduce power bill drastically.

### Standard Accessories Include

Water Level Indicator.  
Stainless Steel Basket.  
Vacuum Breaker.  
Automatic pressure control switch.  
Automatic low water cut off switch.  
Fully automatic PID purging system.