Autoclave Vertical Double Walled - Radial Locking Modular Cabinet Type

Sku: GI1016



www.geomatrixlabs.com

Description

In steam pressure sterilization autoclaves are ideal for all applications requiring, assured and reliable total destruction of all living microorganisms

Our top loading Autoclaves are the ideal companion for the modern microbiologists and institutes, pharma companies, hospitals etc.

FEATURES:

- · In Outer Rectanguar body, inner, outer chamber and support ring are made of Stainless Steel (SS-304 grade)
- · Heaters cover stand and carrier/bucket are made of stainless steel.
- · Lid made of thick machined stainless steel plate.
- · Radial locking lid ensures user safety and safe locking.
- · Lid is fitted with pressure gauge, spring loaded safety valve, double safety valve, manual exhaust valve, vacuum breaker as a safety device and water level indicator.
- · Easy portability & placement, unit is mounted on castor wheels for easy movement
- · Drain valve facilitates easy cleaning of the chamber.
- · Moulded joint less silicon rubber gasket.
- \cdot All internal joints are argon arc welded, ground and polished to give crevice free internals.
- · Miniature circuit breaker for mains on/off plus safety of control panel.
- · Industrial grade energy efficient heaters reduce power bill drastically.

Standard Accessories Include

Water Level Indicator. Stainless Steel Basket.

Vacuum Breaker.

Automatic pressure control switch.

Automatic low water cut off switch.

Fully automatic PID purging system.