

No image found or type unknown

Autoclave Vertical Double Walled – Wing Nut Locking

Sku : GI1008



www.geomatrixlabs.com

Description

GEOMATRIX Vertical Top loading double walled Autoclaves models for high-pressure sterilization are ideal for all applications requiring, assured and reliable total destruction of all living micro-organisms. An Ideal equipment for Laboratories, hospitals, nursing homes, R & D labs, etc. It is used to sterilize instruments by subjecting them to high pressure saturated steam at 121° C for stipulated time.