

No image found or type unknown

Burette Clamp Double

Sku : GE0078



www.geomatrixlabs.com

Description

- Aluminium die-casted body , blue painted with spring loaded
Feature black rubber sleeves for firm gripping of burettes
- Can hold one or two burettes up to 100 ml.
- This one-handed, clamp having scissor-like mechanism provides quick and easy adjustment.
- Insert burette and gently release to grip . Numbers and graduation on burette remain easy to read ,Suitable for retort stands with support rods up to 12.5mm diameter