## Centrifuge Machine Medium – High Speed 10000 R.P.M.

**Sku:** GI1066



www.geomatrixlabs.com

## **Description**

## **Features**

- · Bench Top High Speed Centrifuge Machine
- · Microprocessor based with digital timer 0 -59 minutes.
- · Heavy duty C.R.C. construction finished with powder coating
- · Metalic centrifugal bowl for extra protection and easy cleaning.
- · Fitted with well balanced universal motor (with carbons)
- · 2 lines of 16 character LCD display of r.p.m., Set time and run time
- · Programmable Speed controller from 2000 to 10000 in variation of 100 rpm
- · Last test parameter recall and displays automatically on LCD
- · Lid lock for safe drive
- · Rotors made of Aluminium Alloy to withstand high speed of 10000 rpm

## **TECHNICAL PARAMETER**

Max. Speed: 10000 r.p.m.

Motor: Universal Motor (With carbon)
Power: AC 220V 50 Hz + 10%
Dimensions: 510 x 510 x 390mm
Weight of Machine: 21 kg approx.
Dimension & wt. with Box: 23.5 kg