

No image found for product type unknown

Centrifuge Machine Medium – High Speed 10000 R.P.M.

Sku : GI1066



www.geomatrixlabs.com

Description

Features

- Bench Top High Speed Centrifuge Machine
- Microprocessor based with digital timer 0 -59 minutes.
- Heavy duty C.R.C. construction finished with powder coating
- Metallic centrifugal bowl for extra protection and easy cleaning.
- Fitted with well balanced universal motor (with carbons)
- 2 lines of 16 character LCD display of r.p.m., Set time and run time
- Programmable Speed controller from 2000 to 10000 in variation of 100 rpm
- Last test parameter recall and displays automatically on LCD
- Lid lock for safe drive
- Rotors made of Aluminium Alloy to withstand high speed of 10000 rpm

TECHNICAL PARAMETER

Max. Speed : 10000 r.p.m.
Motor : Universal Motor (With carbon)
Power : AC 220V 50 Hz + 10%
Dimensions : 510 x 510 x 390mm
Weight of Machine : 21 kg approx.
Dimension & wt. with Box : 23.5 kg