

No image found type unknown

## Kipps Apparatus

**Sku :** GE6165



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Description

Geomatrix Labs Presenting Kipp's Apparatus, constructed from polyethylene and comprising three main sections. The top portion features a thistle funnel and a bulb structure, facilitating the pouring of dilute acid to react with a metal sulphide at the base. The evolved H<sub>2</sub>S gas collects in the middle dome and can be released through a leakproof stopcock fitted with a self-lubricating PTFE plug. The base includes a residual outlet fitted with a cap for effortless cleaning.