

No image found or type unknown

MCT Twin Rack (Rack for 0.5 ml & 1.5 ml MCTS)

Sku : GE6251



www.geomatrixlabs.com

Description

Introducing MCT boxes, expertly molded from polycarbonate, ensuring both strength and autoclavable performance. One size can accommodate 80 MCTs of 1.5 ml, while the other holds 100 MCTs of 0.5 ml. These boxes are designed to withstand temperatures ranging from -190°C to +121°C. The transparent cover allows a clear view of the contents inside and is numbered for efficient inventory control.