

No image found for product type unknown

## Meker Burner

**Sku :** GE0114

**Metalware Type 1 0 Col1 :** GE0114A



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Product Image Gallery



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Description

Cast iron heavy base epoxy coated finish, chrome plated brass tube top dia 25mm with a perforated mesh steel head plate which provides small multiple hot flames.

Meker-type burners are in usage when an extremely hot flame is required.