

No image found for product type unknown

## Natural Rubber Tubing Geo Star

**Sku :** GE2003



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Product Image Gallery



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Description

- Made from natural rubber (N.R). In addition to softness, Geo-Excel Grade II tubing has more strength and resistant to dilute acids than Grade III
- The reduced wall thickness of this tubing makes it more flexible and economical for use in Laboratory. Shore hardness 45A, Temperature range -20°C to +70°C, Tolerance  $\pm 0.2$