

No image found for type unknown

PCR tube rack with Hinges – Autoclavable

Sku : GE6258



www.geomatrixlabs.com

Description

Removable hinged lid fits both sides Flip this polypropylene reversible rack to the size of tubes you need. PCR side of the rack holds 8- or 12- tube strips or individual tubes. The opposite side of the rack holds 0.5 mL tubes. Both sides of the rack have 12 wells for 1.5 mL tubes. Autoclavable. Rack measures 206x 131 x 54 mm. Pack of 5.