

No image found for product type unknown

Pipette Pump

Sku : GE6021



www.geomatrixlabs.com

Description

Geomatrix Labs introduces a rapid release pipetting device, designed for accurate pipetting and quick liquid dispensing. The knurled thumbwheel allows effortless drawing and dispensing of liquids, while the releaser enables rapid emptying. The molded chuck, featuring a silicon rubber collar, securely holds pipettes of standard sizes. It offers excellent resistance to acids and alkalis and can be easily disassembled for thorough cleaning.