

No image found for product type unknown


## Platinum Cured Silicon Geo Prime Platinum

**Sku :** GE2002L



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Product Image Gallery



### Description

Platinum Cured Silicon tubing suitable for use in the laboratory. Suitable temperature range  $-50^{\circ}\text{C}$  to  $+250^{\circ}\text{C}$ . Shore hardness 45A. Resistance to strong alkalis, Autoclavable. Excellent chemical and solvent resistance. Very Flexible and tear resistant. Transparent tubing made from non-toxic ingredients.