

No image found

## Powder Funnel

**Sku :** GE6181



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Description

Geomatrix Powder Funnels fabricated from Polypropylene material, featuring a parallel stem to minimize powder bridging and external ribbing for preventing air locks. This shatter-proof alternative to glass offers excellent chemical resistance to various lab chemicals and reagents, making it a valuable addition to lab safety programs. Moreover, these funnels are steam autoclavable at 121°C for efficient sterilization.