

No image found type unknown

## Rubber Stopper Two Hole

**SKU :** GE2025Z



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Product Image Gallery



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Description

Shore hardness 45-50A, Temperature Range -30°C to 70+°C. Based on Natural Rubber (NR). Extrasoft Grade II Light weight & soft, it is resistant to dilute acid, alkalies, ammonia Usage in laboratory.