


No image found

Test Tube Stand (Wire Pattern) 13 mm x 90 Tubes

Sku : GE6398



www.geomatrixlabs.com

Description

Test Tube Stand (Wire Pattern) economical, submersible, and autoclavable racks, a perfect alternative to wire racks. Stackable when empty, they come in various colors and are made from a special blend of polypropylene, ensuring stability in agitated water. These space-saving racks can be easily disassembled for cleaning. The base features a grid with square openings, along with fold-up sides and end bars that securely lock into the rack tops.