

No image found type unknown

## Tissue Culture Plate with Cover

**SKU :** GE6192



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Description

For easier handling, our cell culture plates have rounded edges. A lid with a low profile and less evaporation ensures even development. Select from non-treated or tissue-treated surfaces. The surface of an untreated growth medium is comparable to the surface of a bacteriological Petri dish and more hydrophobic than tissue culture medium. Each plate has tiny feet on the bottom corners, which prevent the microscope lens from being scraped.