

No image found type unknown

Universal Reagent Reservoir 50 ml

SKU : GE6438



www.geomatrixlabs.com

Description

Universal Reagent Reservoir offers three variations: a graduated single reagent reservoir with a maximum volume of 50 ml, a reservoir with 8 separate numbered 5 ml reagent channels, and another reservoir with 12 separate numbered 5 ml reagent channels. Crafted from polypropylene, these reagent reservoirs are autoclavable and compatible with all multi-channel pipettes. They provide convenient and precise dispensing for various laboratory applications.