

No image found type unknown

Water Distillation Unit Double Stage All Quartz

Sku : GE3011



www.geomatrixlabs.com

Description

- Quartz Water Distiller has been redesigned for better performance keeping in view customer need. All Quartz Double Distillation consists of demountable upper and lower boiler
- The lower boiler has a built in heater enclosed in a quartz boiler with a cup on top
- The upper boiler has a built in quartz coil and condensing unit which is mounted on the boiler with the help of ball and socket joint.
- Double walled condenser ensures separate condensation of vapors coming from the two boilers using a single cooling circuit.