

No image found type unknown

## Water Distillation Unit Single Stage

**SKU :** GE3001



[www.geomatrixlabs.com](http://www.geomatrixlabs.com)

### Description

This water distillation unit consists of a boiler made from high purity quartz with built in heater and bottom discharge joint for easy cleaning of deposits.

Spiral condenser made from borosilicate glass 3.3 is fitted on the boiler with receiving adapter.

All unit is mounted on a sturdy powder coated stand.

This unit comes with Safety Control Unit (Optional) to protect the Glass parts in case of water failure or overheating.