

Silicon Rubber Tubing – Geo Prime, I.D 15 mm, O.D 23 mm, Cat No GE2001W



Description :

Silicone rubber tubing suitable for use in the laboratory. Resistance to strong alkalis, Autoclavable.

Excellent chemical and solvent Resistance. Transparent tubing made from non-toxic ingredient. Very

Flexible and tear resistant

Permissible temperature range -50°C to +250°C,

Shore Hardness 45 A

Silicone Rubber Tubing Size I.D 15 mm, O.D 23 mm.