No image found Image not found or type unknown

Rack For Micro Centrifuge (Folding) – 100 places For 1.5 ml MCT’s

Sku : GE6404



Description

Even when fully loaded, racks can be stacked and have special label slots on either side. There are 100 spaces for 1.5 ml MCTs.

These two-tier, multipurpose polypropylene racks can be used dry, in water baths, or in freezers. They are autoclavable. Alpha-numeric indices are given on these racks for simple sample identification. Even when fully loaded, racks can be stacked and have special label slots on either side. There are 100 spaces for 1.5 ml MCTs.